

DIN: 00763136

NET CONTENTS / CONTENU NET: 400 g

Blood Stopper

VETERINARY USE ONLY

USAGE VÉTÉRINAIRE SEULEMENT

ACTIVE INGREDIENTS / INGREDIENT ACTIFS:

Ferrous Sulfate:	840 mg / g	Sulphate de Ferrous
Tannic Acid:	10 mg / g	Acide Tannique
Ammonium Alum:	50 mg /g	Ammonium Alum

Manufactured By:

Equivet Ltd.
61 Longsdale Dr.
London, ON
N6G 1T4

Distributed By:

Can-Vet Animal Health Supplies Ltd.
4-61 Malcolm Rd.
Guelph, ON
N1K 1A7

INDICATIONS: Blood Stopper is used to stop bleeding of minor cuts and wounds.

DOSAGE & ADMINISTRATION: Apply liberally to the wound. After bleeding has stopped, apply again.

CAUTION: KEEP OUT OF REACH OF CHILDREN. For external use only. Keep the container tightly closed. Tetanus antitoxin should be given as a precaution.

INDICATIONS: Blood Stopper est utilisé pour arrêter les coupures et blessures bénignes de saigner.

MODE D'EMPLOI: Appliqué libéralement sur l'endroit affecté. Après le saignement arrêté appliqué une deuxième fois.

CAUTION: TENIR HORS DE LA PORTÉE DES ENFANTS. Pour usage externe seulement. Tenir la flacon bien fermé. Donner ne injection antitétanique par mesure de précaution.